



ogi[®]

Marcamos la diferencia

OGI T-STREET

The most practical transfer wheelchair

Check out our website: www.orthoglobalideas.com

Follow us in our **social media**.



Ultra-light Transfer wheelchair

TECHNICAL DATA

Chassis	Aluminium	Rear wheels	PU 11 inch
Net weight	10,8 kg	Total width	55 cm
Max. user weight	100 kg	Total height	91 cm
Seat height	46 cm	Total depth	91 cm
Seat width	43 cm	Folded size	55 x 26 x 72 cm
Seat depth	42 cm	Packaging size	59,5 x 32 x 75,5 cm
Front wheels	PU 6 inch		

Fast & easy
folding system



Compact design



Self-blocking
handle



Inward folding
footrest



270° foldable
armrest



Removable
& washable
upholstery



10,8 kg

OGI T-STREET

Extendable handle model



Easy transport
and storage

Transport
bag included



ogi.[®]

Marcamos la diferencia



ortho global ideas
Marcamos la diferencia

Ortho Global Ideas
C/ Mas la Mora, 13
08500 España

www.orthoglobalideas.com
info@orthoglobalideas.com
tlf. +34 938 86 26 89