

OGI SMART TILT FUNCTION V-2

Foldable wheelchair
with powered seat tilt
and reclining backrest.



ogi.

Marcamos la diferencia

OGI WHEELCHAIRS

Manual & electric models

Check out our website: www.orthoglobalideas.com

Follow us in our **social networks.**



OGI CHIKA 400



OGI CHIKA 400

FEATURES

- 7 inch front wheels
- 16 inch rear wheels
- Double brake: for user and carer
- Foldable armrest and footplate
- Split backrest

TECHNICAL DATA

ARMREST HEIGHT	74 cm ± 5 mm
BACKREST HEIGHT	45 cm ± 5 mm
SEAT HEIGHT	48 cm ± 5 mm
SEAT DEPTH	44 cm ± 5 mm
SEAT WIDTH	43 cm / 45 cm ± 5 mm
TOTAL WIDTH	58 cm ± 5 mm
TOTAL LENGTH	91 cm ± 5 mm
TOTAL HEIGHT	90 cm ± 5 mm
PACKAGING SIZE	82 x 28 x 78 cm
FOLDED SIZE	80 x 27 x 77 cm
TOTAL WEIGHT	11,4 kg
MAX USER WEIGHT	131 kg



Scan the QR code
FOR MORE INFORMATION

OGI CHIKA 550



OGI CHIKA 550

FEATURES


- 7 inch front wheels
- 22 inch rear wheels
- Double brake: for user and carer
- Foldable armrest and footplate
- Split backrest

TECHNICAL DATA

ARMREST HEIGHT	74 cm ± 5 mm
BACKREST HEIGHT	46 cm ± 5 mm
SEAT HEIGHT	49 cm ± 5 mm
SEAT DEPTH	44 cm ± 5 mm
SEAT WIDTH	43 cm / 45 cm ± 5 mm
TOTAL WIDTH	65 cm ± 5 mm
TOTAL LENGTH	100 cm ± 5 mm
TOTAL HEIGHT	92 cm ± 5 mm
PACKAGING SIZE	91 x 31 x 80 cm
FOLDED SIZE	90 x 30 x 79 cm
TOTAL WEIGHT	12 kg
MAX USER WEIGHT	131 kg



Scan the QR code
FOR MORE INFORMATION



OGI FUN 400



OGI FUN 600QR



OGI FUN 400 FEATURES

- 7 inch front wheels
- 16 inch rear wheels
- Double brake: for user and carer
- Foldable armrest and footplate
- Foldable and removable footrest
- Double cross brace
- Split backrest

TECHNICAL DATA

ARMREST HEIGHT	73 cm ± 5 mm		
BACKREST HEIGHT	48 cm ± 5 mm		
SEAT HEIGHT	49 cm ± 5 mm		
SEAT DEPTH	42 cm ± 5 mm		
SEAT WIDTH	40-42 cm	43-45 cm	46-48 cm
MAX USER WEIGHT	120 kg	131 kg	150 kg
TOTAL WIDTH	59 cm	62 cm	65 cm
TOTAL LENGTH	94 cm ± 5 mm		
TOTAL HEIGHT	92 cm ± 5 mm		
PACKAGING SIZE	72 x 30 x 79 cm		
FOOTRESTS	1,6 kg		
TOTAL WEIGHT*	12,8 kg		
TRANSPORT WEIGHT**	11,4 kg		

* Including footrests, soft cushion and brake handles
** Chassis without footrests.

OGI FUN 600QR FEATURES

- 7 inch front wheels
- 24 inch rear wheels
- Double brake: for user and carer
- Foldable armrest and footplate
- Foldable and removable footrest
- Quick-Release wheel system
- Double cross brace
- Split backrest

TECHNICAL DATA

ARMREST HEIGHT	74 cm ± 5 mm		
BACKREST HEIGHT	48 cm ± 5 mm		
SEAT HEIGHT	49 cm ± 5 mm		
SEAT DEPTH	42 cm ± 5 mm		
SEAT WIDTH	40-42 cm	43-45 cm	46-48 cm
MAX USER WEIGHT	120 kg	131 kg	150 kg
TOTAL WIDTH	62 cm	65 cm	68 cm
TOTAL LENGTH	105 cm ± 5 mm		
TOTAL HEIGHT	92 cm ± 5 mm		
PACKAGING SIZE	78 x 34 x 81 cm		
FOOTRESTS	1,6 kg		
QR WHEELS	3,6 kg		
TOTAL WEIGHT*	13,6 kg		
TRANSPORT WEIGHT**	8,4 kg		

* Including QR wheels, footrests, soft cushion and brake handles
** Chassis without footrests or wheels.



Scan the QR code
FOR MORE INFORMATION



Scan the QR code
FOR MORE INFORMATION

OGI URBAN I 400



OGI URBAN II 300



OGI URBAN I 400 FEATURES

- 7 inch front wheels
- 16 inch rear wheels
- ≤6 km/h max speed
- Range ≥25 km depending on user weight and terrain.
- Foldable armrest
- Split backrest

TECHNICAL DATA

NET WEIGHT	23 kg		
MAX USER WEIGHT	110 kg		
ARMREST HEIGHT	21 cm ± 5 mm		
BACKREST HEIGHT	43 cm ± 5 mm		
SEAT HEIGHT	49 cm ± 5 mm		
SEAT DEPTH	40 cm ± 5 mm		
SEAT WIDTH	41 cm	43 cm	46 cm
TOTAL WIDTH	56 cm	58 cm	60 cm
TOTAL LENGTH	88 cm ± 5 mm		
TOTAL HEIGHT	93 cm ± 5 mm		
BATTERY	Ion Lítio 1 * 28 Ah		
MOTOR	2 * 150W		
PACKAGING SIZE	80 cm * 39 cm * 76 cm		

OGI URBAN II 300 FEATURES

- 8 inch front wheels
- 12 inch rear wheels
- ≤6 km/h max speed
- Range ≥18 km depending on user weight and terrain.
- Foldable armrest
- Removable footrests with tilting and height adjustable footplates
- Split backrest

TECHNICAL DATA

NET WEIGHT	33 kg		
MAX USER WEIGHT	100 kg		
ARMREST HEIGHT	74 cm ± 5 mm		
BACKREST HEIGHT	43 cm ± 5 mm		
SEAT HEIGHT	50 cm ± 5 mm		
SEAT DEPTH	43 cm ± 5 mm		
SEAT WIDTH	41 cm	43 cm	46 cm
TOTAL WIDTH	57 cm	60 cm	63 cm
TOTAL LENGTH	103 cm ± 5 mm		
TOTAL HEIGHT	93 cm ± 5 mm		
BATTERY	Ion Lítio 2 * 11,6 Ah		
MOTOR	2 * 200W		
PACKAGING SIZE	78 cm x 33 cm x 93 cm		



Scan the QR code
FOR MORE INFORMATION



Scan the QR code
FOR MORE INFORMATION

PRODUCT

FEATURES

Aluminium chassis

Foldable design for easy transport

Net weight 42 kg

Max. user weight 150 kg

Lithium batteries 24W/2*6Ah (up to 3 uds.)

Unfolded dimensions: **104 x 59 x 127 cm**

Up to 30 km range

as per number of batteries (variable depending on terrain and user weight)

Electromagnetic brakes

Adjustable and removable headrest

Handles slopes up to 6°

2 brushless motors 250W

Max. speed 8 km/h (depending on terrain)

Folded dimensions: **76 x 59 x 62 cm**

8 inch front wheels

12 inch rear wheels

Equipped with **anti tippers**

COMPACT FOLDING



UP TO 3 BATTERIES



SOFT-CUSHION



ALUMINUM CHASSIS



POWERED
**RECLINING
BACKREST**
90-130°



POWERED
SEAT TILT
0-20°



ortho global ideas
Marcamos la diferencia

Ortho Global Ideas
C/ Mas la Mora, 13
08500 España

www.orthoglobalideas.com
info@orthoglobalideas.com
tlf. +34 938 86 26 89